



# How to Print Sample Collected by / Date & Time in Lab Reports

Please follow the following steps for Sample By

## 1) Setup the tuning in master to show Sample By, Please see the attached



## Step 2: Now create the Sample By Master



### Step 3: Now this Sample By will be in list while accept test

Lab Master - [Test Booking]

Masters Booking Test Perform Printing Transaction Lab Reports Settings Utilities Window Help

Patient ID: 0014129 Date: 05 / 10 / 2015 \* Time (hh:mm): 14 : 3 Lab No: vice.35149  
 Name: Mr. 15s1s Mobile No.:  
 Sex: MALE Date Of Birth: 05 / 10 / 2015 **Sample By: ABC**  
 Referred By: Panel Code:  
 Panel ID:  
 e-mail:

Test ID	Test Name	Rate	Tax(%)	Tax Amt
185	PLATELET COUNT	500.00	00.00	00.00

Total Tests: 1  
 Test Amt.: 500.00  
 Concession: 00.00  
 Home Collection: 00.00  
 Tax Amt: 00.00  
 Net Amt: 500.00  
 Balance: 00.00

Home Collection: 0 Pay Type: CASH Paid: 500 Receipt No.: R-00013681

F1: Select Test F2: Save Entry F3: Perform Test F4: New Booking F5: Payment Type F6: Delete Entry F7: Conc. F8: Add Row  
 F9: Tests Selected F11: Print Receipt F12: Patient Details Ctrl+F7: Payment Mode Ctrl+F2: Show Concession \* Time in 24 Hours  
 Format (Click on Time Label for more Information) Value in this field is coming from the Masters.  
 Ctrl+F3: Clinical Details Alt+F5: Visible Fields Shift + F12 : Print Labels

Start 3 Wind... APMRD... APLab - ... AP - Mic... Untitled... Custome... Lab Ma... untitled... STEP2.J... 2:10 PM

### Step 4: For adding the Sample Collection Time, in the Perform Test press on Additional Information and setup values there

Lab Master

Masters Booking Test Perform Printing Transaction Lab Reports Settings Utilities Window Help

Welcome

**Performing Test**  
 TEST NAME: HAEMATOL  
 PLATELET COUNT

**Additional Informations**  
 Image No.:  
 Ref. By:  
 Performed By:  
 athologyREP. TIME:  
 Sample Collection Date / Time \* 21 / 11 / 2011 10: 50  
 Mandral Reference No.:  
 Type Of Sample:  
 Source of Blood:  
 Investigation:  
 Remarks F2: Ok Esc: Cancel

\* Put time in 24-Hours format. Put Hours and Minutes as 00, if you dont want to print date only

F5: Calculator F7: Print F8: Preview

Accurate Info Soft Pvt.

Start 3 Wind... APMRD... APLab - ... AP - Mic... Untitled... Custome... Lab Master STEP3.J... STEP2.J... 2:11 PM

## Now the Lab Reports will appear like this

Test Report Of vice.35149

100% 1 of 1 F2: Print Email PDF Esc: Exit

Preview

LAB No.	vice.35149	Date	05/Oct/2015
Name	<b>Mr. S15s1s</b>	Sex	M
Sample By	ABC	Sample Date	21/Nov/2011 10:50

**HAEMATOLOGY ROUTINE**

Test Name	Result	Units	Reference Range
PLATELET COUNT	<b>1</b>	10 <sup>3</sup> /mm	(150-450)

----- End of Report -----