

How to repair Access Database of Labmaster


We can repair Access database by two ways.

1. Using Microsoft access
2. Using ODBC Data source

Here I am using you one by one repair method


1. How we can repair with Microsoft Access

We have to select the database Labmaster from running drive



We can open the above file using double click, it will opened in Microsoft Access

In Microsoft 2003 and Microsoft 2007 version we can repair from ---



Using this option


Once you run this option system will repair automatically, after finishing this step close Microsoft Access

And run your labmaster software



- 2 .using ODBC Data source

IN Windows XP Start →Setting →Control Panel → Administrative Tools →ODBC Data source

Microsoft Vista, Microsoft Windows 7 and Window 8


Start →Control Panel\All Control Panel Items\Administrative Tools → 

This screen will be opened after clicking on the Add Button




Double Click on Highlighted Option

Driver do Microsoft Access (*.mdb) → Double click



Click on Repair




Select the path of labmaster folder Using Drives option

Double click on Labmaster Folder or Specified folder there labmaster Database exist

On left hand side you can see the database just double click or press OK

Message will occur after repairing the database



Now close all application

Now your database is repaired and ready to use